

Acellular Dermal Matrix

MegaDerm[®]

www.lncbio.co.kr

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

Freeze/dried type products: e-beam 멸균 실시
Hydration type products: Gamma 멸균 실시

Indication of MegaDerm

Life and Creativity

L&C BIO

Expansion of business field through development of new medical devices

MegaDerm® HD

- Pre-hydrated human acellular dermal matrix

Sub-pectoral breast reconstruction with MegaDerm® HD

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MegaDerm® STAR

- Pre-formed design for pre-pectoral breast reconstruction
- Pre-hydrated human acellular dermal matrix

Pre-pectoral breast reconstruction with MegaDerm® STAR

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For sub-pectoral breast reconstruction

Size (cm)	Thickness (mm)
5 x 14	1.0 - 1.5
6 x 14	
6 x 16	
7 x 18	
8 x 16	

Size (cm)	Thickness (mm)
5 x 16	1.5 - 2.3
6 x 14	
6 x 16	
7 x 18	
8 x 16	

For pre-pectoral breast reconstruction

Size (cm)	Thickness (mm)
12 x 16	1.0 - 1.5 or 1.5 - 2.3
10 x 20	
14 x 14	
16 x 16	
18 x 18	

For pre-pectoral breast reconstruction (MegaDerm Star)

Size (cm)	Thickness (mm)
16 x 16	1.0 - 1.5
18 x 18	
20 x 20	

Volume replacement with diced MegaDerm® in breast conserving surgery

- Ref) Volume replacement with diced acellular dermal matrix in oncoplastic breastconserving surgery : a prospective single center experience, Gwak et al. World Journal of Surgical Oncology, (2020) 18:60

The excision cavity was then filled with the diced MegaDerm® pieces.


Diced type

Size (cm)	Thickness (mm)
5 x 6	5
3 x 6	

- + A, B) Mammography showed diced-MegaDerm as a mass with a well circumscribed margin on the left mediolateral oblique view and left cranio-caudal view.
- + (C) The diced-MegaDerm was observed to have the same echogenicity as fibroglandular parenchyma on B-mode ultrasonography.
- + (D, E) The diced-MegaDerm showed low signal intensity on T2-weighted imaging and no enhancement on contrast-enhanced magnetic resonance imaging.

Satisfaction with cosmetic and overall surgical outcomes in patients

The responses in the cosmetic and overall satisfaction questionnaire revealed that more than 90% of the patients were strongly satisfied.

- + Patients : 120 patients aged 20–80 years with breast cancer who desired BCS between December 2017 and August 2018
- + Results : 1) Cosmetic outcome ($p = 0.250$)
Patient group: 9.7 (± 0.8)
Surgeon group: 9.7 (± 0.8)
2) Overall satisfaction ($p = 0.001$)
Patient group: 9.4 (± 1.0)
Surgeon group: 9.5 (± 1.1)
- + Analysis : Patient and surgeon satisfaction with the surgery outcomes 6 months postoperatively

Score	N (%)	Satisfaction Level										Average score
		Strongly satisfied	Satisfied					Neutral			Dissatisfied	
		10	9	8	7	6	5	4	3	2	1	
Cosmetic outcome	Surgeon	91(77.8)	16(13.7)	8(6.8)	1(0.9)	0	0	1(0.9)	0	0	0	9.7
	Patient	92(78.6)	18(15.4)	5(4.3)	0	1(0.9)	1(0.9)	0	0	0	0	9.7
Overall outcome	Surgeon	78(66.7)	23(19.7)	9(7.7)	3(2.6)	2(1.7)	1(0.9)	1(0.9)	0	0	0	9.4
	Patient	83(70.9)	23(19.7)	6(5.1)	0	4(3.4)	1(0.9)	0	0	0	0	9.5

Postoperative complication rates in breast cancer patients after volume replacement with diced MegaDerm®

Postoperative complication rates in 117 breast cancer patients after volume replacement 6 months postoperatively

- + Reoperation rate 8.5%, Removal case: 2case
- + Hematoma 5.1%, Seroma 6.0%, Fat necrosis 3.4%

Score	Seroma	Red breast syndrome	Infection	Hematoma	Wound edge necrosis	Fat necrosis	Total
Incidence	7(6.0%)	3(2.5%)	3(2.5%)	6(5.1%)	1(0.9%)	4(3.4%)	24(20.5%)
Reoperation	0	1(0.9%)	0	4(3.4%)	1(0.9%)	4(3.4%)	10(8.5%)
ADM Removal	0	1(0.9%)	0	1(0.9%)	1(0.9%)	0	2(1.7%)



MegaDerm®: Cross-linked human acellular dermal matrix

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● Angiogenesis: 6 months after implantation

● Angiogenesis: 6 months after implantation

* Ref) Lee, Ju Hee, Hyung Goo Kim, and Won Jai Lee. "Characterization and tissue incorporation of cross-linked human acellular dermal matrix." Biomaterials 44 : 195-205.

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Products of MegaDerm in Rhinoplasty


Camouflaging of silicon

Dorsum (block type)

Dorsum (carving type)

Tip

Type	Size (cm)	Thickness (mm)
Camouflaging of silicon	1 x 5	0.7 - 0.9
Dorsum (block type)	1 x 5	1.8 - 2.3
		2.3 - 3.0
		3.0 - 4.0
		4.0 - 5.0
Dorsum (carving type)	1 x 5	3.0 (W3)
		4.0 (W4)
		5.0 (W5)
Tip	1 x 1	1.5 - 2.3
		2.5 - 3.0
		3.0 - 3.5
		3.5 - 4.0
		4.5 - 5.0
		5.0 - 6.0

Clinical Study

Comparative skin evaluation after split-thickness skin grafts using 2 different acellular dermal matrices to cover composite forearm defects.


* Lee, Myung Chul, et al., The Journal of hand surgery 42.4 : 297-e1.

Functional Skin Assessment in MegaDerm

AU, arbitrary unit. The MegaDerm grafted area and adjacent normal skin had comparable elasticity and transepidermal water loss but the MegaDerm area was less moist, darker, and redder (mean SD). *P< .05

Variable	MegaDerm Graft Area	Control Area	P Value
Elasticity	0.60 ± 0.15	0.61 ± 0.16	.22
Humidification	24.64 ± 13.48	39.90 ± 14.78	<.05*
Transepidermal water loss, g/h/m ²	10.93 ± 6.38	12.89 ± 6.22	.15
Melanin value (AU)	228.00 ± 86.50	154.00 ± 59.15	<.05*
Erythema value (AU)	324.00 ± 91.00	259.00 ± 81.06	<.05*

Comparative skin assessment results using the ratio of measurements

Aesthetic Evaluation Based on Patient's Point of View

Aesthetic Grade	MegaDerm®	STSG only
Excellent	8	3
Good	15	6
Fair	6	12
Poor	0	2
Total	29	23


Clinical Study

Preventive Effect of Human Acellular Dermal Matrix on Post-thyroidectomy Scars and Adhesions: A Randomized, Double-Blinded, Controlled Trial.

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


Erythema index(AUs)


Clinical Study

Effect of human acellular dermal matrix (MegaDerm®) on infraauricular depressed deformities, Frey's syndrome, and first bite syndrome following parotidectomy: a multi-center prospective study.

* Joo Hyun Kim, et al. Gland Surg 2021;10(2):670-677.

The implantation of MegaDerm following total parotidectomy.

- + A. Surgical field after total parotidectomy with preservation of the facial nerve.
- + B. MegaDerm was designed and inserted between flap in parotid bed and residual parotid tissue.

Evaluation of Frey's syndrome of parotidectomy patients in the MegaDerm-use group and control group

	Month 3		Month 6		Month 12	
	ADM	Control	ADM	Control	ADM	Control
Patients (%)	3 (4.9%)	14(25.9%)	4 (6.5%)	12(22.2%)	4 (6.5%)	8 (14.8%)
Total score (0-10)	0.23±1.2	3.48±1.54	0.36±1.31	3.57±1.1	4 (6.5%)	3.64±0.85
Frey Qx (1-4)	1.05 ± 0.34	2.13±0.48	1.25±0.24	2.36±0.78	1.32±1.15	1.41±0.78
P-value	0.032 / 0.027 / 0.095		0.037 / 0.024 / 0.315		0.018 / 0.024 / 0.287	

Visual analogue scale for the contouring deformity and subjective satisfaction score of parotidectomy patients in the MegaDerm-used group and control group.

Countouring deformity	Month 3		Month 6		Month 12	
	MegaDerm	Control	MegaDerm	Control	MegaDerm	Control
Not detectable deformity	31	8	25	5	22	4
Mild deformity	26	23	32	22	34	23
Moderate deformity	4	15	4	16	5	16
Severe deformity	0	7	0	10	0	10
Extremely severe deformity	0	1	0	1	0	1
Subjective satisfaction score	4.26±0.51	3.78±0.97	4.16±0.32	3.63±0.92	3.96±0.54	3.57±0.86
P-value (VAS/SSS)	0.001 / 0.001		0.001 / 0		0.001 / 0	

Clinical Study

Surgical Outcomes of Tympanoplasty Using Acellular Dermal Allograft: A Randomized Controlled Study.


* Kim, Sung Huhn et al. ACTA OTORHINOLARYNGOLOGICA ITALICA 2018;38:554-562;

Compared to autologous graft materials, MegaDerm is an effective alternative as a TM graft material with similar graft success rates and postoperative hearing results

* TM: tympanic membrane

Representative tympanic membrane


- + A. Perforation of the tympanic membrane(TM)
- + B. The use of MegaDerm as a TM graft
- + C. Otoscopic views of the TM at 6 months postoperatively with complete TM closure


Hearing results: air-conduction, bone-conduction, air-bone gap

- + Fig A, B) Preoperative and postoperative hearing results of the perichondrium group and MegaDerm group measured at 6 months (A: perichondrium group, B: MegaDerm group) postoperatively. The air-conduction thresholds significantly improved in both groups, whereas the bone-conduction thresholds remained stable.
- + Fig C) Changes in air-bone gaps (ABGs) after tympanoplasty at 6 months postoperatively. Significant reductions in ABGs were observed in both groups, and the amount of reduction of ABGs and postoperative ABGs were not significantly different between the two groups.


Perichondrium group

MegaDerm group

Perichondrium & MegaDerm

Replacement of pericranium: Wound depression protectant after craniotomy

Clinical case:
protectant of depression after mini-craniotomy

- + A) Frontotemporal Depression :
Pterional Craniotomy without MegaDerm
- + B) No Frontotemporal Depression :
Pterional Craniotomy with MegaDerm

Size (cm)	Thickness (mm)
4 x 5	0.6 - 0.8
	1.0 - 1.5
6 x 8	0.6 - 0.8
	1.0 - 1.5


Arthroscopic rotator cuff repair with MegaDerm®

Clinical Benefits


- + Augmentation of soft tissue repairs
- + Biologically versatile scaffold for ready incorporation to support regeneration
- + Protects delicate tissue with a durable and safe matrix

Case	Size (cm)	Thickness (mm)
Massive tear case : Augmentation	4x5, 3x4, 3x3	2.0
Ratear case	4x5, 3x4, 3x3	2.0
Small tear case	1x3	0.7 - 0.9


1. Massive tear case : Augmentation

2. Retear case


3. small tear case: interposition


Pterional mini-Craniotomy

Intraoperative photograph after pterional mini-craniotomy


Intraoperative photograph after MegaDerm® insertion

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Ref) Characterization and tissue incorporation of cross-linked human acellular dermal matrix, Biomaterials 44 : 195-205.

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번호	논문 제목	저널명	일자
1	Lee, Ju Hee, Hyung Goo Kim, and Won Jai Lee. "Characterization and tissue incorporation of cross-linked human acellular dermal matrix."	Biomaterials 44 : 195-205.	2015 Mar
2	Kang, Sang-Wook, et al. "Preventive Effect of Human Acellular Dermal Matrix on Post-thyroidectomy Scars and Adhesions: A Randomized, Double-Blinded, Controlled Trial."	Dermatologic Surgery 41.7: 812-820.	2015 Jul
3	Kim, Jae-Young, et al. "Esthetic neck dissection using an endoscope via retroauricular incision: a report of two cases."	Journal of the Korean Association of Oral and Maxillofacial Surgeons 40.1: 27-31.	2014 Feb
4	Jeong, Woon Hyeok, et al. "Acceleration of osteogenesis by platelet-rich plasma with acellular dermal matrix in a calvarial defect model."	Child's Nervous System 32.9: 1653-1659.	2016 Sep
5	Park, Gui-Yong, et al. "Acellular dermal matrix as a core strut for projection in nipple reconstruction: approaches for three different methods of breast reconstruction."	Archives of plastic surgery 43.5 :424.	2016 Sep
6	Suh, Man Koon, et al. "Augmentation rhinoplasty with silicone implant covered with acellular dermal matrix in Asian noses"	Journal of Craniofacial Surgery 28.2 : 445-448.	2017 Mar
7	Lee, Jun Ho, et al. "The effect of sterile acellular dermal matrix use on complication rates in implant-based immediate breast reconstructions."	Archives of plastic surgery 43.6 : 523.	2016 Nov
8	Lee, Myung Chul, et al. "Comparative skin evaluation after split-thickness skin grafts using 2 different acellular dermal matrices to cover composite forearm defects."	The Journal of hand surgery 42.4 : 297-e1.	2017 Apr
9	Kim, Jiye, et al. "Use of Acellular Allogenic Dermal Matrix (MegaDerm) in Orbital Wall Reconstruction: A Comparison With Absorbable Mesh Plate and Porous Polyethylene."	Journal of Craniofacial Surgery 28.7 : e644-e649	2017 Oct
10	Kim, June-Kyu, and Yong Seong Kim. "Comparison of the Amount of Drainage in 3 Acellular Dermal Matrices in Implant-Based Breast Reconstruction: A Retrospective Study."	Archives of Aesthetic Plastic Surgery 23.2 : 68-72.	2017

번호	논문 제목	저널명	일자
11	Yang, Chae Eun, et al. "Usefulness of Cross-Linked Human Acellular Dermal Matrix as an Implant for Dorsal Augmentation in Rhinoplasty."	Aesthetic plastic surgery 42.1 : 288-294.	2017 Nov
12	Kim, Sung Huhn et al. "Surgical Outcomes of Tympanoplasty Using Acellular Dermal Allograft: A Randomized Controlled Study"	ACTA OTORHINOLARYNGOLOGICA ITALICA 2018;38:554-562;	2018 Dec
13	Kim, Chang-Hoon, and Sang Chul Park. "Homologous Tissue for Dorsal Augmentation."	Facial Plastic Surgery Clinics of North America 26(3):311-321	2018 Aug
14	Seo H. Lee et al. "Dual Coverage of the Inferior Pole with Conjoined Fascial Flap and Acellular Dermal Matrix for Immediate One-Stage Breast Reconstruction with a Prosthetic Implant"	Aesthetic Plast Surg. 2018 Oct;42(5):1213-1219.	2018 Apr
15	Lee, D., Kim, Y. S., Roh, T. S., & Yun, I. S. (2019). Cryptotia recurrence lowering technique with additional acellular dermal matrix graft. Archives of craniofacial surgery, 20(3), 170.	Arch Craniofac Surg Vol.20 No.3, 170-175	2019 Jun
16	A comparative study between sterile freeze-dried and sterile pre-hydrated acellular dermal matrix in tissue expander/implant breast reconstruction	3.63±0.92	2019 May
17	Volume replacement with diced acellular dermal matrix in oncoplastic breastconserving surgery: a prospective singlecenterexperience	World Journal of Surgical Oncology	2020 Mar
18	Comparison of Volume Retention and Biocompatibility of Acellular Dermal Matrix/Hyaluronic Acid Filler to Autologous Fat Grafts in a Mouse Model	Aesthetic Plastic Surgery	2020 Mar
19	The relationship of human acellular dermal matrix thickness on complication rate and patient-reported outcomes in implant based immediate breast reconstruction	Gland Surg 2021;10(1):90-100	2021 Jan
20	Effect of human acellular dermal matrix (Megaderm™) on infra-auricular depressed deformities, Frey's syndrome, and first bite syndrome following parotidectomy: a multi-center prospective study	Gland Surg 2021;10(2):670-677	2021 Feb

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