



BT-150
Hi-bebe
Super

CONTENTS

- INTRODUCTION
- FEATURES & BENEFITS
- WHAT'S INCLUDED
- SPECIFICATIONS



Hi-bebe Super



Intended use of Breast pumps



- Breast pump are intended to be used by lactating women for expressing breast milk in home environments.
- Expressing milk by breast pump is indicated to alleviate the symptoms of engorgement.
- In case of symptoms of mastitis, breast pumps may support the healing process by removing breastmilk from the affected breast.

Hi-bebe super usage

- To relieve sore and cracked nipples
- Separated from baby because of work or travel
- To bring out flat or retracted nipples
- When baby cannot suck mother's nipple because of severe nipple retraction
- To store the breastmilk because the milk expression amount is irregular



Hi-bebe Super Benefits & Features

Various Levels for Expression & Massage Mode

BPA Free/ Latex Free Materials -> Good to baby and mother

Nursing Night Lamp -> can be used during night.

Mobile Application -> Compatatible with Andoroid & Apple


Closed pumping system -> prevention of breast milk backflow



What's Included



Carrying bag



Mobile Application



Bottle Sets



Main Unit



Funnel Sets(21/24/27/32 mm)

Specifications

- Vacuum: 20~290 mmHG
- Weight: 1.7kg
- Warranty: 1 Year
- Size(Main body): 190(W) X 190(D) X 79(H) mm
- Rechargeable battery: Operating Time: 120min
Charging Time: 150min
- Operating environment: Temperature 5 ~40 C
Relative humidity: 0 ~95% non-condensing
Atmospheric pressure: 80 ~106 kpa
- Options : Carrying bag, Car charging cable, Funnel, Nursing Night lamp, Mobile App
- Certificates: KFDA, CE, FDA, ISO13485, GMP



THANK YOU



Tel : 82-31-750-0340

Fax : 82-31-750-0344

E-mail : bistos@bistos.co.kr

[Http://www.bistos.co.kr](http://www.bistos.co.kr)

© Copyright 2001 Bistos Co., Ltd. All Rights Reserved. This material is protected by Korea and other copyrights and may not be copied, sold or redistributed in any form without the written permission of Bistos Co., Ltd.